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Dear Customer,

thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please read through these **instructions for use** carefully. If you have any questions, please contact your doctor or medical retailer.

Intended purpose

SecuTec Lumbo is a medical device. It is an orthosis for lordosis reduction and relief of the lumbar spine with a bridging frame.

Indications

- Most severe lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Most severe (pseudo) radicular lumbar syndrome
- Spondylolysis / spondylolisthesis, grades I–III
- Moderate to severe lumbar deformity in cases of facet syndrome / osteoarthritis
- Lumbar spinal canal stenosis
- Moderate to severe degeneration / muscular insufficiency
- Vertebral fractures (lumbar spine) with vertebral body damage
- Condition following herniation, e.g. post-operative
- Condition following discectomy

Risks of using this product

Caution*

- SecuTec Lumbo must only be worn in accordance with the specifications contained in these instructions and used for the areas of application (places of use) listed.
- SecuTec Lumbo is a prescribed product that should be applied under a physician's guidance.
- To ensure the best possible fit, the SecuTec Lumbo must be individually adapted by a trained specialist¹ to optimally suit the patient's requirements (adaptation of the bridging frame).
- Only a trained specialist may carry out the (initial) fitting of the orthosis and provide instruction.
- Use in combination with other products, e.g. as part of compression treatment (compression stockings), must be discussed beforehand with your physician.
- The SecuTec Lumbo is to be worn under your T-shirt or undershirt.
- Should you notice any unusual changes (e.g. increased symptoms), please contact your physician without delay.
- Incorrect use of the support will cancel any product liability.
- Do not allow the product to come into contact with ointments, lotions or substances containing grease or acids.
- Side effects involving the body as a whole have not been reported to date. Correct fitting is assumed. Do not put the SecuTec Lumbo on too tightly as this could lead to local pressure symptoms. In rare cases, constriction of blood vessels and nerves may occur.

 Contains natural rubber latex (elastodiene), which can cause skin irritation in persons with an allergy. Take off the orthose if you notice any allergic reaction.

Contraindications

Hypersensitivity reactions harmful to health have not been reported to date. In the following conditions such aids should only be put on and applied after consultation with your physician:

- Skin disorders / injuries in the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness and excessive heat build-up
- Impaired sensation and movement disorders, e.g. in the case of diabetes mellitus
- Severe restriction of the performance of the heart and lungs (risk of a rise in blood pressure if the aid is worn during intense physical exertion)

Application instructions

- A –Knitted fabric portion
- B –Fastening with finger pocket
- C –Tensioning straps

- D –Pocket with bridging frame
- E –Abdominal pad

Fitting

The SecuTec Lumbo must be adapted by a trained specialist¹ before the first fitting.

1 Position the orthosis centrally, at waist height, so that the lower back is covered. To ensure optimal effectiveness, the orthosis should be positioned as low as possible without pressing uncomfortably against the groin area while the patient is sitting upright. **2** Then slide your fingers from the sides into the finger pockets provided on the fastening and steadily pull the fastening forward.

3 Then press the left half of the fastening onto your stomach and position the right half of the fastening on top of the left half so that it can be fastened to it on the left-hand side. When fastening the two halves, pull your fingers out of the left-hand finger pocket first and then out of the right-hand one. **4** Then take hold of the two tensioning straps and pull them forward, steadily and at the same time, until the desired level of tension (stabilization) is achieved.

5 Then fasten the two tensioning straps close together on the large fastening. **6** + **7** Now fasten the two upper tensioning straps.

Removal

To take off the orthosis, first unfasten the four tensioning straps. Loosely refasten the straps onto the two halves of the fastening, then open the large fastening. Remove the orthosis, then refasten the large fastening to avoid damage.

Cleaning instructions

Note: Never expose the SecuTec Lumbo to direct heat (e.g. heater, sunlight, in your car etc.!) This may cause damage to the material, which can impair the effectiveness of the SecuTec Lumbo.

The sewn-in label of the the SecuTec Lumbo, containing information about the product name, size, manufacturer, washing instructions and CE identification, is located on the inside of the knitted material.

Before washing, remove the bridging frame from the velour pocket, unfasten the abdominal pad and close all Velcro fastenings to avoid damage.

Wipe the bridging frame clean with a mild detergent. Wash the knitted fabric portion and the abdominal pad separately by hand with a mild detergent at 30 °C.

Allow the bridging frame, abdominal pad and knitted section to air dry.

When reinserting the bridging frame into the velour pocket provided, make sure that the two markings ("up" and "outward") are correctly oriented.

Part of the body this product is used for

In accordance with the indications (lower back). See Purpose.

Maintenance instructions

With correct handling and care, the product is practically maintenance-free.

Assembly and fitting instructions

The SecuTec Lumbo is supplied in a standard size and ready for use. The bridging frame must be individually adapted by a trained specialist and positioned in the knitted section.

Technical specifications / parameters

The SecuTec Lumbo is a orthosis for the lower back. It consists of a knitted section, a velour pocket for inserting the bridging frame, an abdominal pad and six double-sided micro Velcro fastenings. The velour pocket and bridging frame are available as spare parts on request.

Information about re-use

The SecuTec Lumbo is intended for your treatment only. It has been adapted specifically for you. For this reason, do not pass the SecuTec Lumbo on to other people.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty.

The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the SecuTec Lumbo have not been observed.

Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the specialist were not observed
- The product was modified arbitrarily

Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Polyamide (PA), Polyurethane (PUR), Stainless steel, Viscose (CV), Glass fiber (GF), Elastodiene (ED), Polyethylene (PE), Polyoxymethylene (POM), Polyether foam, Ethylene propylene diene monomer rubber (EPDM), Polyester (PES), Carbon fiber (CF), Cotton (CO), Polypropylene (PP), Elastane (EL), Copolyester (TPC-ES)

 – Medical Device

 – DataMatrix Unique Device Identifier UDI

Specialists¹

General notes

Check the optimal and individual fit of the SecuTec Lumbo after the first fitting. Adapt the bridging frame (CFR frame) according to the area of application (indications) and the anatomy of the patient. Practice the correct procedure for putting on the orthosis with the patient.

The product must not be changed in any improper way. Failure to comply with this requirement may adversely affect product performance, thereby voiding any product liability.

The SecuTec Lumbo is intended for the treatment of one patient only.

Part of the body this product is used for

In accordance with the indications (lower back). See Purpose.

Set-up and assembly instructions

- A –Knitted fabric portion
- B –Fastening with finger pocket
- C –Tensioning straps
- D –Velour pocket for bridging frame (CFR frame)
- E –Abdominal pad
- F –Bridging frame (CFR frame)
- G –Double-sided micro Velcro fastenings

Caution*

Wear protective gloves to avoid burns when handling the hot-air device.

To prevent damage to the material, please pay attention to the following: Use round bending hooks only. Do not heat the CFR frame at temperatures higher than 190 °C. The optimum temperature for manipulating the frame is between 170 °C and 190 °C. Maintain a minimum distance between the hot-air device and the CFR frame. Heat the CFR frame evenly and over the entire surface (not just at certain points).

Bend the CFR frame in large arches only (no sharp bends). You can immerse the CFR frame in cold water to speed up the cooling process. This will not have any negative impact on the properties of the material.

You can heat and reshape the material more than once. The heating process does not visually alter the surface of the material. In order to find out the appropriate time to bend the material, test it during the heating process.

1 Heat the CFR frame with a hot-air device and shape it to produce the necessary bridging angle in accordance with the area of application (indication) and to suit the patient's anatomy. Make sure that the "up" marking on the frame is correctly oriented. Check the fit of the orthosis on the patient in both a standing and a sitting position after each modification.

2 Lay out the velour pocket in front of you. The sewn-in label indicates which side should face upward ("up"). Then turn the velour pocket so that the opening is pointing toward you. **3** Place the CFR frame in the pocket and close it. **4** Place the knitted section in front of you with the sewn-in label on the top right. Secure each pair of double-sided micro Velcro fastenings onto the velour strip on the knitted section. The fastenings must be spaced out at distances that enable the velour pocket to be securely connected to the frame and the knitted section.

5 Fasten the velour pocket to the bridging frame and the knitted section.

6 Fasten the abdominal pad onto the right-hand half of the fastening in front of you.

Practice the process of putting on and taking off the orthosis with the patient.

Version: 2021-03

* Note on the risk of personal injury (risk of injury, harm to health and accidents) or damage to property (damage to the product).

¹ A specialist is any person who is authorized according to the state regulations for fitting and instruction in the use of orthoses which are relevant to you.

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