BAUERFEIND AG
Triebeser Straße 16
07937 Zeulenroda-Triebes
Germany
P +49 (0) 36628 66-40 00
F +49 (0) 36628 66-44 99
E info@bauerfeind.com
BAUERFEIND.COM

MOS-Genu

Montageanleitung
Assembly instructions
Instructions de montage
Instrucciones de montaje
Istruzioni per il montaggio

MOS Genu
103 648

de	deutsch	2	sr	srpski	20
en	english	3	sl	slovenski	22
fr	français	4	ro	romanian	23
nl	nederlands	5	el	ελληνικά	24
it	italiano	7	tr	türkçe	25
es	español	8	ru	русский	26
pt	português	9	et	eesti	28
sv	svenska	10	lv	latviešu	29
no	norsk	11	lt	lietuvių	30
fi	suomi	12	uk	українська	31
da	dansk	14	kk	қазақ	32
pl	polski	15	ja	日本語	33
cs	česky	16	ko	한국어	34
sk	slovensky	17	ar	عربي	35
hu	magyar	18	zh	中文	36
hr	hrvatski	19			

de deutsch

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Only for specialists¹

Dear retail partner,

these **assembly instructions** contain important information on the adjustment options of the MOS-Genu.

In addition, please read and observe the **instructions for use** carefully. If you have any questions, please contact your Bauerfeind sales partner.

Procedure for First Time Application

1. Preparation

- Loosen the 4 adjustment cams above and below the knee ring.
- Mark the theoretical anatomical knee joint centre at 30° flexion.
- Please choose the size of the knee joint padding for the individual. MOS-Genu comes ready-fitted with post-op pads on the inner side of the joint for post-operative use.

2. Static fitting

⚠ IMPORTANT!

The shell sections and the straps must lie flat on the skin. Do not use sharp-edged tools to contour the orthosis. Most adjustments can be performed manually.

- The initial fitting is carried out without the knee pads.
- The MOS-Genu is placed against the knee joint in 30 degree flexion such that the mechanical pivot point of the MOS-Genu corresponds with Nietert's compromise axis.
- Check that the support shells lie flat against the knee both vertically and circumferentially; adjust if necessary.
- A round adjusting hook with safety cap should be used for major adjustments.
- Select the knee pads of the required thickness, fasten and turn the Velcro-backed pads until flat contact with the knee is achieved. The post-op pads must be removed before the knee padding is attached using the Velcro straps.
- Centre the distal (lower) popliteal fossa strap with the integrated knee pad, tighten and secure on both sides.
- Take the proximal (upper) popliteal fossa strap and proceed as for the distal strap.
- Now don, center and use the Velcro to attach the wide straps distally and proximally (from below and above). For the long version, finally attach the straps of the extension shells using Velcro.
- The brace is held in the individual bowleg / knock-knee position whilst the straps are attached. The adjusting cams should now be screwed tight to maintain this position. **Important!** Do not attempt to move adjusting cams once tightened, otherwise the material will be stressed.
- On the long version, the height of the distal (lower) and proximal (upper) extension shells can be adjusted to the full extent of the lever arm.
- If necessary, the protective shinbone padding, also supplied, can be affixed in order to relieve pressure pain in the shinbone area - or if swelling is present in the area of the tibial tubercle.

3. Adjustment examples

- Two stop pins can be set to limit joint flexion and extension.
- The extent of permitted movement can be read off, in degrees, in the window on the outer aspect of the knee.
- Flexion / extension can be restricted by placing the stop pins at the appropriate positions.
- Angles of 0°, 20°, 40°, 60° etc. can be set with the stop pin in the long slot and angles of 10°, 30°, 50° etc. with the stop pin in the short slot.
- The orthosis is shipped with a 0°/0°/90° flexion / extension initial setting.


Fig. 1: Schematic diagram of possible joint settings

4. Dynamic fitting

- Check the effect of MOS-Genu on the patient's movement.
- For gait anomalies, loosen the adjustment cams and in 20 steps let the orthosis run in, then retighten adjustment cams. Important! Do not attempt to move adjusting cams once tightened, otherwise the material will be stressed.
- Repeat this procedure until the desired gait is reached.
- TORX screws have been used exclusively on the MOS-Genu. A matching screwdriver is also supplied. TORX screws allow for high tightening torques as it is virtually impossible for the screwdriver to slip out of the TORX drive.

5. Bowleg / knock-knee settings

- The graduated adjusting cams allow the brace to be placed in the precise bowleg / knock-knee position of the leg.
- Loosen the cam screws and set the guide cams first (identified by spots in the illustration). The follower cams will turn in the opposite direction.
- Precise adjustment is achieved by alternating the rotation of both cams distally (lower) and proximally (upper).
- The distal (bottom) and proximal (top) scale settings on the cam pairs should be the same. Adding the distal (bottom) and proximal (top) readings together will give the individual setting.
- In appropriate applications, targeted relief of the ligaments on the inner or outer aspects of the knee can be achieved by correcting the bowleg / knock-knee setting.

Covering disk	Upper (proximal)	Lower (distal)	Setting
Bandy legs	3°	4°	△ = 7°

Fig. 2: Varus setting of 7° overall on the left leg

6. Dismantling the long version

- Depending on the application and the progress made by the patient in the rehabilitation phase, it is possible to simply dismantle the brace and turn it into the short version.
- Release and remove the screws on the distal (lower) and proximal (upper) support shells.
- Pull the distal (lower) and proximal (upper) extension shells out of the support shells.
- Compensate for the resulting hollows on the inner aspects of the shells with the padded strips provided.
- Place the orthosis in position and check that the fit is correct.

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- A 4-point system is implemented with the aid of the two popliteal fossa pads (Fig. 3), e.g. for an anterior instability (drawer sign).
- For a posterior instability (drawer sign) or back-knee (genu recurvatum), there is a long popliteal cross-strap available to replace the short popliteal straps.
- Twist the cross-strap through 180° and thread through the pad loops. Pull through the fixing loops on the distal (lower) joint ring and cross over in a proximal (upper) direction and secure (Fig. 4).
- The popliteal cross-strap can also be used after surgery to make the brace more comfortable to wear.


Fig. 3: ACL Treatment

Fig. 4: PCL Treatment

MD – Medical Device

UDI – DataMatrix Unique Device Identifier UDI

Version: 2026-04

1 A professional is any person who is authorized according to the state or institutional regulations as well as other provisions relevant to them for fitting and instruction in the use of orthopedic aids. As the manufacturer, we strongly recommend that this person be trained in orthopedics or possess comparable acquired competencies as well as dispose of orthopedic infrastructure.

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