AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

Bauerfeind Benelux B.V.
Waarderveldweg 1
2031 BK Haarlem
THE NETHERLANDS
P +31 (0) 23 531-9427
F +31 (0) 23 532-1970
E info@bauerfeind.nl
BELGIUM
P +32 (0) 2 527-4060
F +32 (0) 2 792-5345
E info@bauerfeind.be

BOSNIA AND HERZEGOVINA

Bauerfeind d.o.o.
Meše Selmovića 19
71000 Sarajevo
P +387 (0) 33 710-100
F +387 (0) 33 619-422
E info@bauerfeind.ba

CHINA

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 english

Dear patient,

Thank you for choosing a Bauerfeind product. We wish you a speedy recovery. Please read and observe these instructions for use closely and in their entirety. If there is something you do not understand or if you have any questions, do not use the product until you have consulted your physician or professional¹ again. Keep these instructions for use in case you want to look something up later.

Intended purpose

OmoTrain is a medical device. It is a support for the relief and stabilization of the shoulder joint.

Indications

- Chronic, post-traumatic or post-operative irritation
- Prevention / recurrence prevention
- Osteoarthritis (joint wear)

Contraindications

For the following clinical pictures, usage of the product is only indicated upon consultation with your physician:

- Skin disorders/lesions on the part of the body to be treated, especially if inflammatory manifestations as well as raised scars with swelling, redness and overheating are present
- Impaired sensation and movement disorders of the arms (e.g. diabetes mellitus)
- Impaired lymphatic drainage – including unspecified soft tissue swelling located in distal direction of the fitted support
- Venous occlusion in the arm (arm vein thrombosis)

Ask your physician whether your conditions fall into these categories.

Notes before first usage

If you wish to use your Bauerfeind product for acute symptoms, injuries or healed previous conditions, you urgently need to obtain professional medical advice and observe it before using the product. Restrictions may apply to you which you must take into account and which the professional¹ will clarify.

In case of identified or suspected problems with blood pressure and intense physical exercise, prolonged use of the product can have additional effects on your blood pressure. After you have been prescribed/you have purchased the product, only use it according to the indication and under observance of any additional instructions given by your physician/ professional¹. If you wish to use the product in combination with medication, other medical devices or other treatments, consult your physician/ professional¹ beforehand. If any additional conditions or injuries occur during the period of application of the product, also consult your physician immediately.

Risks of application

- Do not arbitrarily modify the product as this may cause it to not help as expected and/or worsen your health.
- The product must not be applied too firmly as this may cause local pressure symptoms or, in rare cases, constrict underlying blood vessels or nerves. If this happens, loosen the strap on the product.
- If you notice any changes or an increase in symptoms related to the application of the product (such as tissue congestion, skin irritation, swelling, etc.), pause the application and consult your physician/ professional¹.
- If you completely doff the product, treatment is insufficient and additional damage to the part of the body affected may occur during additional strain. Be especially careful when this is the case.
- Additionally, restrictions due to the indication may apply which the professional¹ will explain to you. In particular, driving vehicles, other means of transportation or machines may only be possible under restrictions. The same applies to the use of ladders, tools and sports equipment. When in doubt, we recommend to suspend these activities.
- Persons requiring supervision due to their age or mental/physical impairment must be supervised during every application of the product.
- Avoid contact of the product with fatty and acidic agents, ointments, lotions, bases, sand or salt water.
- If hypersensitivity has been observed with you, consult your physician before the first application of the product. If hypersensitivity/ an allergic reaction occurs in relation with the application of the product, pause the application and consult a physician immediately.

- Never expose your product to direct sunlight, heat or cold for longer periods of time.

Contains natural rubber latex (elastodiene) which may cause skin irritations if hypersensitivity is present. Doff the product if you notice any allergic reactions.

Application instructions

- Observe the instructions for donning and doffing of the product as well as any additional instructions for use given by your physician/ professional¹ under any circumstances.
- Since the product unfolds its effects especially during physical activity, doff it during longer rest phases.
- Check the product for integrity before every use. Only intact products may be used.
- Make sure to have your physician/ professional¹ show you how to don/doff the product and explain to you how it works until no doubts remain. If the affected part of the body changes in diameter, the product may have to be readjusted/newly prescribed.

Donning of the product

A – Knitted section of the support, B – Velcro straps, C – Eyelet, D – Message pad (Delta pad)

- 1** Fasten one of the OmoTrain straps and pull it out through the eyelet. To avoid damaging your clothes or the OmoTrain, fasten the ends of the strap loosely onto the strap underneath it. **2** Secure the knitted section onto the fastened straps. The opening should be pointing toward you. **3** Then guide the affected arm through the knitted section of the support until the shoulder is well covered. The message pad (Delta pad) should be resting on the bony protuberance at front and back of the shoulder. If this is not the case, adjust the position of the support or the message pad. **4** Arrange both straps in a flat (without twisting), diagonal position across the stomach. **5** Take hold of them by inserting your hand through the loop from the inside. **6** Pull the loop behind your head and slip your arm through it. **7** Guide the rear strap across your back toward the affected arm and secure it in place. **8** Then guide the front strap around your abdomen toward the affected side. Guide the strap that you unfasted to begin with back through the eyelet and fasten it. **9** Then tighten the abdominal strap until the required degree of stabilization is achieved. Fasten the end of the strap.

Doffing of the product

Undo both straps. To avoid damaging your clothes or the OmoTrain, fasten the ends of the straps loosely. Slip your arm out of the sling first, and then slip it out of the knitted section of the support.

Cleaning instructions

Refasten the message pad after washing. To do this, lay the OmoTrain with the straps facing upward on the long side of the knitted fabric. Fasten the message pad onto the knitted section with the long tapered end pointing toward you. Observe the supplementary instructions on the sewn-in label/sticker on your product. Wash the product at 30 °C/86 °F. Regular care ensures the ideal effectiveness of the product.

Area of use

Shoulder joint

Notes on reuse

The product is intended for the individual treatment of one patient only.

Warranty

The legal provisions of the country in which the product was purchased apply. If you suspect a warranty claim, please first contact the person from whom you purchased the product directly. The product must be cleaned before making a warranty claim. If the instructions for handling and caring for the product have not been observed, the warranty may also be affected or excluded.

The warranty does not cover:

- Off-label use
- Failure to observe the instructions of the physician/professional¹
- Arbitrary product modifications

Disclaimer

In general, do not self-diagnose, self-treat or self-medicate unless you are a professional¹. Before first using our medical device, actively consult a physician or professional¹ since this is the only way to assess the individual effect of the product on your body and determine any possible risks of application/side effects. Please observe the advice of this professional¹ as well as all notes and instructions in this document or its online version (including texts, images, graphics, etc.), including any excerpts.

Obligation to report

Due to regional legal regulations, you are obliged to report any serious incidents involving the use of this medical device to both the manufacturer and the competent authority for medical devices (in Germany: Bundesinstitut für Arzneimittel und Medizinprodukte - Federal Institute for Drugs and Medical Devices).

Disposal

Please dispose of the product in accordance with local regulations.

Material content

Polyamide (PA), Styrene-Ethylene-Butylene-Styrene copolymer (TPE SEBS), Elastane (EL), Elastodiène (ED), Viscose (CV), Polyoxyméthylène (POM), Polyuréthane (PUR), Cotton (CO), Polyester (PES)

 – Medical Device

 – DataMatrix identifier as UDI

Date of information: 2026-04

1 A professional is any person who is authorized to adjust supports and orthoses and instruct others in their use in accordance with the applicable state and institutional regulations and other guidelines. As the manufacturer, we strongly recommend that this person disposes of orthopedic technical training or comparable acquired skills, as well as orthopedic technical infrastructure.

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