

# ORDER FOR CUSTOM-MADE PRODUCTS

## GenuTrain®



BAUERFEIND.COM

141

141

Identifier \_\_\_\_\_  
(e.g. patient name, abbreviation) \_\_\_\_\_

Fax Germany +49 (0) 800 001 05 25  
Fax Austria +43 (0) 800 44 30 131  
Fax Switzerland +41 (0) 5 64 85 82 59

Quotation

Order

	Circumference	Length
	f <sup>1</sup>	e - f <sup>1</sup>
	f	e - f
	e	
	d	e - d
	c	e - c
	b <sup>°</sup>	e - b <sup>°</sup>

**Knee**

Right  Left

**Color**

Titan  Beige

**Extras**

Siliconized edge\* (add 5 cm to the length)

Velcro fastening\*

Top\*\*

**Patient details\*\*\***

Male  Female

Obese  Pronounced muscles

Age \_\_\_\_\_ Height \_\_\_\_\_

\* Subject to additional charges

\*\* Not possible in combination with a siliconized edge; not circular

\*\*\* Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

**All measurements must be entered to ensure proper treatment.**

f<sup>1</sup> Thigh: Measure the circumference at least 14 cm above the middle of the knee (e). If the circumference at point f<sup>1</sup> is more than 47 cm, then the circumference must be measured at a higher point, where e - f<sup>1</sup> is at least 17 cm. (note: length e - f<sup>1</sup> standard GenuTrain = 19 cm)

f Base of thigh (approx. middle of e and f<sup>1</sup>)

e Middle of knee

d Narrowest point below the knee (approx. 5-6 cm below the knee)


c (biggest) Calf circumference

b<sup>°</sup> lower end  
Length e - b<sup>°</sup> at least 12 cm  
(note: length e - b<sup>°</sup> standard GenuTrain = 13 cm)

**Notes:**

- Take measurements with a slight flexion (30°).
- The length is measured on the outside of the leg.  
The distance between e - f<sup>1</sup> should not be less than the distance between e - b<sup>°</sup>.
- If the circumference of f<sup>1</sup> exceeds 60 cm, we recommend attaching a silicone top band.

**The purchaser is responsible for ensuring that all size information is correct. Custom-made solutions cannot be returned.**



Stamp

Customer no. \_\_\_\_\_

Signature \_\_\_\_\_

Date \_\_\_\_\_