


ORDER FOR CUSTOM-MADE PRODUCTS

GenuTrain® P3

BAUERFEIND.COM

Identifier _____
 (e.g. patient name, abbreviation) _____

	Circumference	Length
	f¹	short: e - f¹ = 18 cm long: e - f¹ = 22 cm
	f	e - f
	e	
	d	e - d
	c	e - c
	b¹	short: e - b¹ = 15 cm long: e - b¹ = 17 cm

All measurements must be entered to ensure proper treatment.

- f¹ Thigh
 Short version: Measure the circumference 18 cm above the middle of the knee (e).
 Long version: Measure the circumference 22 cm above the middle of the knee (e).
- f Base of thigh
- e Middle of knee
- d Narrowest point below the knee
- c Calf circumference
 Short version: Measure the circumference at a maximum of 12 cm below the middle of the knee (e).
 Long version: Measure the circumference at a maximum of 14 cm below the middle of the knee (e).
- b¹ Calf circumference
 Short version: Measure the circumference 15 cm below the middle of the knee (e).
 Long version: Measure the circumference 17 cm below the middle of the knee (e).

Notes:

- Take measurements with a slight flexion (30°).
- The length is measured on the outside of the leg.
- The selection of the short/long version depends on the individual measurements of the patient. The long version is recommended for patients from a height of 1.70 m.
- Please always specify the position of the corrective guide:
 - Standard (as with standard size range): The corrective strap is located at the lateral (outer) edge of the patella.
 - Alternative option: The corrective strap is located at the medial (inner) edge of the patella.
- The following minimum circumference measurements are required for supplying the product:
 - Applies to short and long versions**
 - f¹ (18 cm or 22 cm above the knee): at least 39 cm
 - f (9 cm or 13 cm above the knee): at least 37 cm
 - e (middle of knee): at least 30 cm
 - b¹ (15 cm or 17 cm below the knee): at least 28 cm

Quotation

Order

Knee

- Right Left
- Short
 (length of support, thigh [e-f¹] = 18 cm
 length of support, lower leg [e-b¹] = 15 cm)
- Long
 (length of support, thigh [e-f¹] = 22 cm
 length of support, lower leg [e-b¹] = 17 cm)

Color

- Titan Nature

Extras

- Siliconized edge* (add 5 cm to the length)

Position of the corrective guide

- Standard Alternative option

Patient details**

- Male Female
 Obese Pronounced muscles

Age _____ Height _____

* Subject to additional charges – please consider the additional length when measuring patients

**Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

The purchaser is responsible for ensuring that all size information is correct. Custom-made solutions cannot be returned.

Stamp

Customer no. _____

Signature _____

Date _____