

ORDER FOR CUSTOM-MADE PRODUCTS


GenuTrain®

BAUERFEIND.COM

Identifier _____
 (e.g. patient name, abbreviation) _____

Quotation

Order

	Circumference	Length
	f ¹	e - f ¹
	f	e - f
	e	
	d	e - d
	c	e - c

Knee

Right Left

Color

Titan Nature

Extras

Siliconized edge* (add 5 cm to the length)

Tightening straps*

Velcro fastening*

Top** Bottom

Patient details***

Male Female

Obese Pronounced muscles

Age _____ Height _____

* Subject to additional charges

** Not possible in combination with a siliconized edge; not circular

***Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

All measurements must be entered to ensure proper treatment.

f¹ Thigh: Measure the circumference at least 14 cm above the middle of the knee (e). If the circumference at point f¹ is more than 47 cm, then the circumference must be measured at a higher point, where e - f¹ is at least 17 cm.

f Base of thigh

e Middle of knee

d Narrowest point below the knee

c Calf circumference: Measure the circumference at least 12 cm below the middle of the knee (e).

Notes:

- Take measurements with a slight flexion (30°).
- The length is measured on the outside of the leg. The distance between e - f¹ should not be less than the distance between e - c.
- If the circumference of f¹ exceeds 60 cm, we recommend attaching a silicone dotted band.



Stamp

Customer no. _____

Signature _____

Date _____