


ORDER FOR CUSTOM-MADE PRODUCTS

MalleoTrain®

BAUERFEIND.COM

Identifier _____
 (e.g. patient name, abbreviation) _____



Quotation

Order

Foot

Right

Left

Color

Titan

Nature

Extras

Velcro fastening*

Measuring points for pads**

① _____

② _____

Patient details***

Male

Female

Obese

Pronounced muscles

Age _____

Height _____

* Subject to additional charges

** Special position of the pads (position of ankle bones differ substantially from the standard measurements)

***Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

All measurements must be entered to ensure proper treatment.

b¹ Base of calf: Measure the circumference at least 16 cm above the heel (A).

b Smallest circumference above the ankle bone

h Instep (circular measurement of the instep over the heel)

a Circumference behind the MTP joint of the little toe

A Heel

A-a Foot length without toe

Notes:

- Only specify the pad measurement points if special positioning is required.
- Distance ① from the medial malleolus to the lateral malleolus (measured horizontally across the Achilles tendon)
- Distance ② from the medial malleolus to the lateral malleolus (measured over the sole of the foot)

Stamp _____

Customer no. _____

Signature _____

Date _____